
# Intra-SWING®

Composite Biface  
*Dualsided Composite*

**IS180 / VENTRO-S**

WWW

Chirurgie Digestive  
*Digestive Surgery*

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MAG x 30

### FACE VISCÉRALE VISCERAL SIDE

Enduction de polyuréthane (PUR) macroperforé  
*Macro-perforated polyurethane (PUR) coating*

Poids Weight - 80 g/m<sup>2</sup>

Porosité Porosity\* - 1.4% (C)

La technologie exclusive PUR

*Exclusive PUR technology*

- Limite les adhérences aux tissus.  
*Reduced risk of adhesions*
- La barrière anti-adhérente permanente offre une protection durable aux viscères.  
*Permanent non-adherent barrier offering a durable visceral protection*

La macroperforation Ø 1 mm

*Ø 1 mm macro-perforation*

- Facilite le drainage des fluides.  
*Ease the fluid drainage.*

Fil de traction en polypropylène monofilament USP#1

*Traction threads in monofilament polypropylene USP#1*

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## CLINICAL AND TECHNICAL DATA

Pre-clinical data:

- PUR : stable, biocompatible and adhesion-free biomaterial1
  - Intraperitoneal placement study (18 rabbits)2
  - Implant extraction at 4, 9 and 13 months : all intact, well integrated and no local sepsis
  - Non-occlusive and loose adhesions in 3/18 rabbits
- INTRA-SWING® prosthesis are clearly harmless and adhesion-free.

Clinical data: Clinical data:

- Intermediary analysis of the INTRA-SWING® observational study (ETOBSI : 71 patients, post-operative follow-up at 2 years at 69%, at 1 year 24%, 7% out of the study) : recurrence 0%, infection 2,8% (2/71).
- Observational clinical study Olmi et al., 20103 (19 patients, post-operative follow-up at 20 months): recurrence 0%, infection 0%
- Clinical market follow-up (117 patients, post-operative follow-up at 1 month): recurrence 0%, infection 0,8%

## Dualsided Composite

Intra-Swing® composite, a complete solution for ventral hernia repair.

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
### CLINICAL AND TECHNICAL DATA

HAS Report, November 20084 :

Clinical data analysis issued from surgeries performed with composite meshes, intraperitoneal placements (medium post-operative follow-up > 2 years)

- e-PTFE visceral side : recurrence ≤ 10%, infection ≤ 9,5%
- Collagen visceral side : recurrence from 2 to 15%
- Two e-PTFE visceral sides : recurrence ≤ 7%, infection ≤ 3,8%

Technical and data :

Issued from internal and external tests (independent accredited laboratories). Our implants respect the recommendations from HAS, European Pharmacopoeia and standards : NF S94-801 and ISO 13993.

## IS180

Dispositif médical classe IIb - Class IIb medical device

IS180 1000 R    Ø 10 cm

IS180 1200 R    Ø 12 cm

IS180 1600 R    Ø 16 cm

IS180 1015 S    10 x 15 cm

IS180 1520 S    15 x 20 cm

IS180 2030 S    20 x 30 cm

IS180 2535 S    25 x 35 cm

## VENTRO-S

Dispositif médical classe IIb - Class IIb medical device

IS180 1200 VO    Ø 12 cm

IS180 1015 VO    10 x 15 cm

IS180 1520 VO    15 x 20 cm

IS180 2030 VO    20 x 30 cm

IS180 2535 VO    25 x 35 cm

SWING  
TECHNOLOGIES

a THT bio-science division

THTBIO-SCIENCE.COM

CONCERN-ENERGOMASH

Exclusive Distributor of THT Bio-Science  
"Concern-Energomash" CJSC  
26/8 Azatutyan ave., Yerevan, Armenia  
Tel.: 011 87 87 17  
info@c-e.am | www.c-e.am